

Earth Science Technology Conference - 2005

**New Laser Techniques for Detection of Radicals,
Isotopes, and Reactive Intermediates from Robotic
Aircraft and Conventional Aircraft for the
Aura Satellite Collaborative Science Effort**

J. G. Anderson, E. J. Moyer, F. N. Keutsch, T. F. Hanisco, and G. S. Engel

June 28, 2005

James G. Anderson

Harvard University


12 Oxford Street, Link Bldg.

Cambridge, MA 02138

Presidential Initiatives

- US Climate Change Science Program: Strategic Plan

1. Reduce uncertainty and related system
2. Understand the ser and managed ecos related global chan
3. Explore the uses ar to manage risks an and change



ow the Earth's climate
e future;

bility of different natural
systems to climate and

s of evolving knowledge
ted to climate variability

- Global Earth Observation Initiative: 2003




NRC Decadal Report Requested by NASA and NOAA


PHASE I

Earth Science and Applications from Space:
Urgent Needs and Opportunities to Serve the Nation
Explicitly Identify the Critical Importance of Societal
Objectives
Emphasize the Required Decision Support Structure

AURA Collaborative Science: Proposed *in situ* Airborne Observations, Hypotheses and Deployment Strategy in Support of UARP



June 2003







Four Primary Focus Areas


- Decadal Forecast of UV Dosage Levels
- Linking Nitrate, Sulfate, Heavy Metal, and Organic Source Strengths, on a Global Scale, to Changes in the Optical Properties of the Atmosphere, a.k.a. Short-Wave Forcing
- The Aura Satellite Collaborative Science Endeavor
- Technology Innovation for the UAV - Small Satellite Era

Effect of the Montreal Protocol




Ozone Loss in Northern Hemisphere from Satellite






Ozone Photochemistry

AER Model Time Slice: 47°N, March 1993




Enhanced Bromine:

- ↘ \uparrow ozone depletion due mainly to $\text{BrO} + \text{ClO}$ cycle
- ↘ $\text{BrO} + \text{HO}_2$ cycle becomes significant O_3 sink below 16 km
($\text{BrO} + \text{HO}_2$ does not drive O_3 depletion if $\text{Br}_y^{\text{TROP}}$ is constant over time)




The Coupled Tropospheric Chemistry,
Carbon Exchange, Climate Problem:
A Low Cost, Integrated Technology,
International, Airborne Approach




WHY THE ARCTIC WARMS QUICKLY

There are many reasons. A key one: as ice melts, darker seawater and exposed land reflect less of the sun's energy, making the remaining ice melt even faster



Proposed Altair Payload



Cloud/Particle
LIDAR

Pressure
Temperature
Telemetry/GPS


Tracer Measurements

Methane (DFG/Herriott)
N₂O (QC/Herriott)
Carbon Dioxide (DFG/Herriott)
Carbon Monoxide (QC/Herriott)
Formaldehyde (Fiber/CRDS)
Water (Lyman Alpha)
Total Water (Lyman Alpha)
Water (DFB/Herriott)
Ozone (UV Absorption)
H₂O/HDO (ICOS or Fragmentation LIF)


Radical Measurements

OH (Ti:Saph/LIF)
HO₂ (Ti:Saph/LIF)
NO₂ (Ti:Saph/LIF)
ClO/ClOOCl/ClONO₂ (RF)
BrO/BrONO₂ (RF)
IO (Ti:Saph/LIF)


1) Temperature control



2) Convective influence



IPCC 4AR dq/dt: Difference Between Model Forecasts



Water Isotopologues: Motivation

- Study of dehydration mechanisms (response to forcing)
- Priority 1 measurement for upcoming science missions
- Development need has been identified
 - Instrument for vapor phase is needed
 - Science requirements demand improved sensitivity
- Direct absorption:
 - simultaneous measurement of H_2O , HDO , H_2^{18}O
- **Substantial pathlength enhancement:**
 - **Integrated Cavity Output Spectroscopy (ICOS):**
> 4km pathlength


HDO measurement requirements

	Science needs
Sensitivity <i>in: Integration time or: Spatial scale</i>	~ 50 ppt 3 s 600 m
Precision	3 % per species (20-25 ‰ in δD)
Accuracy	3 % per species
Hysteresis time	< 3 s

Design of the mid-IR ICOS Flight Instrument for Measurement of HDO, H₂¹⁸O, and H₂O

- **Optical System:**
 - Maintain sensitivity in transfer to robust and compact flight design
- **Gas Sampling System:**
 - Particle free sampling, minimization of wall effects and trapped volumes, constant P in ICOS cell
- **Thermal Management:**
 - Constant T of optical system and of gas volume in cell
- **Electronics:**
 - Compact design for a low pressure environment
- **Software:**
 - Algorithms for robotic control of flight instrument, data acquisition and post-flight data reduction

Integrated Cavity Output Spectroscopy (ICOS)



$$\frac{dl}{dt} = T \times P - l \frac{c(1-R + \alpha \times L)}{L}$$

$$I_{ss} = \frac{PTL}{c(1-R + \alpha \times L)}$$


$$\tau = \frac{L}{c(1-R + \alpha \times L)}$$

$$L_{eff} = \frac{L}{1-R}$$

$$L_{eff} = \frac{L}{1-0.9998} = 5000L$$

- Steady State cw transmission monitored
- Ideal for measuring multiple species
- **Substantial pathlength enhancement: 4200 m**


Off-axis ICOS¹



- On-axis alignment gives high “noise” from cavity resonances
- Off-axis alignment reduces cavity resonances
- Passive cw technique with unsurpassed sensitivity (30 km pathlength)


¹Engel, Keutsch *et al.* in preparation (2005).

Off-axis ICOS




- On-axis alignment gives high “noise” from cavity resonances
- Off-axis alignment reduces cavity resonances
- Passive cw technique with unsurpassed sensitivity (30 km pathlength)

Harvard ICOS instrument



¹Moyer, E.J., Sayres, D.S., Keutsch F.N. *et al.* in preparation (2005).


Flight Modifications of Optical System



- Short optical path and locking optical mounts (minimizes alignment effects of vibrations)
- AR coated lenses and wedged and AR coated QCL window (minimizes fringing effects)
- Removable TE-cooled alignment laser
- Mechanically more stable optical cell
- Spare QCL and spare optics


Motivation: In situ water isotope instruments require validation

- We expect instrument artifacts to contaminate water measurements.
 - Sampling
 - Optical
 - Software analysis
- Independent measurements can help identify artifacts.







H₂O



HDO



‰ D-depletion




Time (s)

Design of ICOS Flight Cell

Flight Modifications: Accommodating Optical and Sampling Considerations

- Novel design of cavity with parallel mirrors and **NO** adjustable or moveable parts (vibration insensitive)
- Thermistors with high impedance and fast time response
- Electropolished and coated surfaces to minimize wall effects
- Pharmaceutical grade connectors and valves, orbital weld joints, P-port placement and cell inlets/outlets designed to eliminate trapped volumes



Inlet Side:


- Rear facing inlet for particle exclusion
- Heaters for gas, feedback to T_{gas} in cell

Outlet Side:


- Large diameter tubing, high throughput
- Scroll Pump: oil-free to protect mirrors
flush time of ca. 3s

ICOS Cell Gas Inlet Design

Inlet ports are located asymmetrically



FLUENT Simulation




Flight Modifications

- Inlet into gas cell results in swirling gas flow
- Pump ports are in center of cell
- Design excludes P gradients and stagnant region in front of mirrors

Fit algorithm yields measurement precision comparable to noise-equivalent absorption

- SNR on the spectrum is > 50 , and both fit residual and derived mixing ratio yield similar values.
- Integration time here is $\frac{1}{2}$ s, 6 x shorter than flight integration times. Projected flight SNR is > 130 in these conditions.
- Increased data acquisition rate and fit improvements should allow further increases in SNR by a factor of 5-10.
- Stratospheric SNRs are more than sufficient for relevant science issues.




Synthetic spectrum in black represents the signal obtained by traditional tunable diode laser absorption spectrometers in the same conditions. The long ICOS optical pathlength produces deeper absorption lines.


Laboratory technique → flight instrument:

- Target molecule: HDO (and H₂O, H₂¹⁸O)
- Platform: NASA's WB57-F high-altitude aircraft
- Test flights: November 2004
- Validation strategy: comparison with simultaneous observations of HDO and H₂O by 3 other independent instruments

HDO chosen because of:

- High scientific importance (priority 1 for Aura science missions)
- No adequate existing technology for stratospheric measurements
- ICOS can meet science needs






Ly- α Water Vapor
 H_2O



HOxotope
 H_2O , HDO



**Total
Payload**



ICOS
 H_2O , HDO, H_2^{18}O , CH_4



Ly- α Total Water
 H_2O (vapor + condensate)



Test flight results

Debut flights for two new water isotope instruments

Independent techniques, complementary strengths

ICOS Isotope Instrument



- Mid-IR absorption spectroscopy using new cavity-based technique
- *(Integrated cavity output spectroscopy)*
- Enhanced sensitivity (x 40) because of 4 km optical path
- Multiple species (H_2O , HDO , H_2^{18}O , H_2^{17}O , CH_4)


HOxotope Instrument



- Photofragment laser-induced fluorescence of OH and OD
- *(Heritage of Harvard HOx instrument)*
- Contamination-free sampling (HO_x radicals lost on wall contact)
- Improved sensitivity (x 2-10) over conventional techniques



H₂O isotopes: Photolysis - Fluorescence detection of HDO/H₂O

Excimer Lamp Photolysis




$$h\nu = 8\text{W}$$

$$\lambda = 172\text{ nm}$$




OH and OD LIF
 $\lambda = 287\text{ nm}$




Radical vs. Molecular Sampling


Molecular water exchanges with walls. OH and OD radicals are lost irreversibly




Consistency between all four instruments




Vertical profiles are plausible




Tropospheric variation is not noise




Model for tropics



Mid-lats observations





HDO in-flight measurement capability

	Science needs	ICOS	HOxotope
Sensitivity <i>in: Integration time or: Spatial scale</i>	~50 ppt 3 s. 600 m	60 ppt 3 s. 600 m	250 ppt 10 s. 2 km
Precision	3 % (20-25 ‰)	3 %	---
Accuracy	3 %	---	---
Hysteresis time	< 3 s.	< 3 s. for cirrus	negligible

Conclusions: ICOS

- Most sensitive airborne mid-IR spectrometer
- Intercomparison with other instruments:
 - Fast time response
 - No evaporation of condensed phase
- Upgrades in progress:
 - Sealed cell
 - Pressure regulation system
- Water isotopologue (gas + condensed phase) measurements allow study of, e.g.:
 - Dehydration, transport, cloud microphysics and response to forcing
- Harvard ICOS instrument ready and available for Aura validation missions

Summary: HOxOTOPE

- Flight data high points:
 - Hoxotope worked as well in flight as in the lab
 - Fast time constant in H₂O and HDO sampling
 - Absence of unexplained artifacts
 - Good relative accuracy
 - Reasonable signal to noise (HDO \pm 250 pptv/4s)
- Post flight schedule:
 - Calibration: 5% (50 per mil)
 - Sensitivity: factor of 5 for moderate effort
 - Investigate total water sampling mode

Development of Miniaturized Intra-Cavity DFG, and Pulsed Fiber-Optic Laser Systems

